

# scRNA-seq analysis pipeline using Snakemake

**MSc Data Science, FHNW**

Crispin M. Lang  
24 November 2025




# Content

- Project goals
- Introduction
- Dataset
- Analysis Pipeline Overview
- Results
- Version control

# Project goals

- Snakemake workflow
  - Use wildcards
  - Add a configuration file
  - Use a mix of shell scripts and Python
  - Create a conda environment
  - Graphical output:
    - UMAP plot
    - DAG
- Version control with GitHub
  - Minimum of 2 commits

# Project overview



# Introduction

## **scRNA sequencing:**

- Transcriptomic technologies
- Exploration of cellular heterogeneity
- high-dimensional dataset.

## **Significance:**

- uncover diversity in heterogeneous cell populations
- grouping cells with similar expression profiles.
- which genes are expressed and in what quantities


# Dataset

## Data:

- pbmc3k (peripheral blood mononuclear cells)
- ~2,700 cells, ~33,000 genes
- different cell types

## Preprocessing:

- min. 200 genes/cell
- min. 3 genes/cell
- 2000 final cells selected
- Final output: AnnData object containing cells (observations) × genes (variables).



# Analysis Pipeline Overview


```
snakemake_project/  
├── snakefile  
├── config.yaml  
├── envs/  
│   └── scanpy.yaml  
├── scripts/  
│   ├── preprocess.py  
│   └── analysis.py  
├── data  
├── results/  
└── README.md
```

# Scanpy:


- streamlines the workflow
- ensures compatibility between steps (same framework)
- AnnData structure keeps data and metadata together
- reduces the complexity (no data moving needed)



# Scanpy - preprocessing



# Scanpy - analysis



# Snakefile:

```
snakemake_project/  
├── snakefile  
├── config.yaml  
├── envs/  
│   └── scanpy.yaml  
├── scripts/  
│   ├── preprocess.py  
│   └── analysis.py  
├── data  
├── results/  
└── README.md
```

```
# Load the configuration file  
configfile: "config.yaml"  
  
# Get list of samples from config  
SAMPLES = config["samples"]
```

# Snakefile:

```
snakefile_project/  
├── snakefile  
├── config.yaml  
├── envs/  
│   └── scanpy.yaml  
├── scripts/  
│   ├── preprocess.py  
│   └── analysis.py  
├── data  
├── results/  
└── README.md
```

```
# config.yaml  
samples: ["pbmc3k"]  
  
# Filtering thresholds for quality control  
min_genes: 200      # filter out cells with fewer than 200 genes expressed  
min_cells: 3        # filter out genes expressed in fewer than 3 cells  
  
# Feature selection  
n_top_genes: 2000   # number of highly variable genes to keep for PCA/UMAP
```

# Snakefile:

```
snakemake_project/  
├── snakefile  
├── config.yaml  
├── envs/  
│   └── scanpy.yaml  
├── scripts/  
│   ├── preprocess.py  
│   └── analysis.py  
├── data  
├── results/  
└── README.md
```

```
rule all:  
    input:  
        # Collect outputs for each sample using list comprehension  
        expand("results/{sample}/umap_{sample}.png", sample=SAMPLES),  
        expand("results/{sample}/adata_{sample}.h5ad", sample=SAMPLES)
```

# Snakefile:

```
snakemake_project/  
├── snakefile  
├── config.yaml  
├── envs/  
│   └── scanpy.yaml  
├── scripts/  
│   ├── preprocess.py  
│   └── analysis.py  
├── data  
├── results/  
└── README.md
```

```
rule preprocess:  
    # Input: path to the 10x Genomics matrix directory for this sample.  
    # We assume the 10x data (matrix.mtx, features.tsv, barcodes.tsv) are in data/{sample}/  
    input:  
        "data/{sample}/"      # directory with 10x data for the sample  
    output:  
        # Save intermediate AnnData after filtering & normalization  
        "results/{sample}/adata_{sample}_filtered.h5ad"  
    params:  
        # Pass filtering parameters from config to the script  
        min_genes=config["min_genes"],  
        min_cells=config["min_cells"]  
    threads: 1  
    conda:  
        "envs/scanpy.yaml"      # Use Scanpy conda environment for this rule  
    script:  
        "scripts/preprocess.py"  # This script will read input and apply preprocessing
```

# Snakefile:

```
snakemake_project/  
├── snakefile  
├── config.yaml  
├── envs/  
│   └── scanpy.yaml  
├── scripts/  
│   ├── preprocess.py  
│   └── analysis.py  
├── data  
├── results/  
└── README.md
```

```
rule analyze:  
    input:  
        # Input is the filtered AnnData from the previous step  
        h5ad="results/{sample}/adata_{sample}_filtered.h5ad"  
    output:  
        # Final outputs: (1) UMAP plot image, (2) final AnnData with all results  
        umap_plot="results/{sample}/umap_{sample}.png",  
        adata_final="results/{sample}/adata_{sample}.h5ad"  
    params:  
        # Pass HVG and other parameters from config  
        n_top_genes=config["n_top_genes"]  
    threads: 1  
    conda:  
        "envs/scanpy.yaml" # Same environment (Scanpy) for this analysis step  
    script:  
        "scripts/analysis.py" # This script performs PCA, clustering, UMAP, etc.
```

```
# envs/scanpy.yaml
channels:
  - conda-forge
  - bioconda
  - defaults


dependencies:
  - python=3.10
  - scanpy=1.9.3
  - anndata=0.9.2
  - numpy<2

  - matplotlib-base
  - python-igraph
  - leidenalg
```

```
snakemake_project/
├── snakefile
├── config.yaml
├── envs/
│   └── scanpy.yaml
├── scripts/
│   ├── preprocess.py
│   └── analysis.py
├── data
├── results/
└── README.md
```




# Running the analysis




# Results

- Identified 6 clusters
- Resolution parameter
- target genes
- Cell type identification
- Explored 2 options




# Results

```
sc.pl.rank_genes_group_heatmap()
```




# Results

`sc.pl.rank_genes_group()`




# Results

`sc.pl.umap()`



# Results



**0: IL32: CD4<sup>+</sup> T cells**

**1: FTL: CD14<sup>+</sup> Monocytes**

**2: NKG7: NK cells**

**3: MS4A1: B cells**

**4: PPBP: Platelets**

**5: LYZ: Dendritic cells**

# GitHub

- Initialized public repository
- Initially committed all data -> failed
- Wrote .gitignore to exclude large files
- Committed after each large change in code structure
- Wrote README for reproducibility
- Final commit today with the pptx slides

# Thank you for your attention!